

# MARSDEN

## HM-201D Automatic Height Measure



- Measuring range: 120cm - 200cm
- Graduations: 1mm
- Height measure for premium column scales.
- Sends data straight to the scale

# Marsden HM-201D: Automatic Height Measure

- Measuring range: 120cm - 200cm
- Graduations: 1mm
- Height measure for Marsden's premium column scale
- Automatically sends height data to the scale
- Simple and easy to use
- Incorporates folding cursor for safety
- Patients height read at eye level
- Covered by 1 year warranty



## Automatically send data

The HM-201D is an automatic height measure which can send data straight to our popular DP-3810 indicator, available on our column scales.

As the indicator on the height measure is lowered to sit atop the patients head, the height data will automatically feed into the scale to help calculate BMI without inputting additional data.

## Pair with a column scale

Our HM-201D is available with our M-100 and M-110 column scales as standard as they feature from our best-selling DP-3810 indicator.

The HM-201D can be attached when the column scale is first purchased, or if you would like to add the HM-201D to your existing height measure, we can guide you through the process.

## Easy to use

Once the height measure is set up, either on the M-100 or the M-110, the height measure is incredibly easy to operate and begin gathering accurate readings.

Simply as your patients to stand on the scale, just as they would to gain their weight reading. Then slowly and carefully move the cursor down until it is sat atop your patients head. The height will then be on screen.

**Measuring Range:** 120cm - 200cm

**Graduations:** 1mm

**Dimensions:** 1000mm x 75mm x 90mm

**Weight:** 0.6kg