

MARSDEN

M-600BT

Marsden Hoist Weighing Attachment with Bluetooth

Lightweight & Portable for use in Home, Community or Hospital

Suits Most Hoist 'Coathanger' Sling Support

Safety Clips to Secure

Bluetooth for Apple and Android devices



MARSDEN



Bariatric Scales & Hoist Weighers

M-600BT

Description

The Marsden M-600BT hoist weighing attachment. This scale has a lightweight design for use on hoists with coat hanger shaped spreader bars, clips on and easy transfer.


The M-600BT is ideal for use within the care home and by the community health visitor. A well-padded and comfortable carry case makes the hoist scale genuinely portable from home to home. The clips on the scale make it very easy to move from hoist to hoist to weigh the patient.

The M-600BT is now MDD approved. They are the only hoist weighing attachments that are legal for use in the UK.

Bluetooth connection to Apple and Android devices.



Features

Capacity	200kg
Graduations	100g
Dimensions	630mm x 275mm x 160mm
Weight of scale	6kg
Power Supply	 Mains / Rechargeable
Adaptor	3pin tip +ve 12V 2A
Bluetooth	V2.1

Reliable. Accurate. Elegant.

