



UK PATENT



US PATENT PENDING



EU PATENT PENDING

# MARSDEN

## M-999

### Patient Transfer Scale



- One device for lateral transfer and weighing of patients
- Weigh immobile patients in seconds: less stress, more dignified
- Administer drug dosages with confidence
- Class III Approved and MDD Approved

[www.marsden-weighing.com](http://www.marsden-weighing.com)



## Weighing critical patients is a crucial process. But it's also a challenging one.

Weighing patients correctly is essential to **ensuring accurate drug dosages** and **administering the correct treatment**.

But for patients requiring urgent, immediate medical attention this is not always easily achievable. **How do you weigh critical and immobile patients accurately, quickly and easily?**

## The Patient Transfer Scale: Welcome to the future of patient weighing

The idea of a transfer board with integrated scale was born of frustration: An acute nurse working for NHS Lanarkshire was determined to find an **easier way to accurately weigh immobile patients**.

The result, the M-999, **revolutionises the weighing of patients in hospitals**. It doesn't just make getting an accurate weight reading easy: it actually makes weighing **part of the already-existing process of transferring a patient from stretcher to bed, or bed to bed**.

This means it's **quick and easy, and saves you time**.

## Get an instant, accurate weight - when you need it most

Healthcare professionals have a duty of care to ensure a safe environment for patients.

When knowing the weight of a patient is absolutely crucial to delivering the right care and treatment, the M-999 **embeds confidence and accuracy in your processes**.

Use the M-999 **throughout the hospital pathway** to enhance outcomes - particularly for trauma patients in emergency departments where **time is crucial**, stroke patients to ensure a quick 'door to needle' time, plus Cardiac, ICU, HTU and X-ray patients.

[www.marsden-weighing.com](http://www.marsden-weighing.com)





## Easy to use means minimal training and effective results

We've made the M-999 **as lightweight as possible** so it's easy to take it to exactly where you need to and get an instant, accurate weight.

**Grab handles** conveniently positioned down each side mean it's easier to move. And because it's all about a simple, precise weight reading, there are no complicated functions to get your head around.

## Key features

- MDD Approved/Class III Approved
- Powered by rechargeable battery
- Large, bright weight display
- Grab handles for portability
- Easy to clean - wipe clean surface
- Hold function



## Tested by 30 hospitals - and we've shaped the design to answer everyone's need

As part of the initial research and development of the M-999, **over 30 hospitals** saw, used and fed back their thoughts on the scale. Trialling the scale as it would be used in practice, and getting feedback from **relevant healthcare professionals and the wider multi-disciplinary team**, contributed massively to the design and functionality of the production version.


Subsequent intensive testing has meant the M-999 delivers precision weighing to Class III medical standards.

## Specification

- Capacity/graduations: 250kg x 500g
- Power: Rechargeable battery
- Dimensions: 1800mm x 700mm x 30mm
- Weighing area dimensions: 1800mm x 510mm

[www.marsden-weighing.com](http://www.marsden-weighing.com)

# Patient Transfer Scale: Key features at-a-glance



## MARSDEN

Tel: +44 (0) 1709 364296

E-mail: [sales@marsdengroup.co.uk](mailto:sales@marsdengroup.co.uk)

[www.marsden-weighing.com](http://www.marsden-weighing.com)  
[www.patienttransferscale.com](http://www.patienttransferscale.com)

**f** MarsdenWeighing

**t** @MarsdenWeighing

**YouTube** [www.youtube.com/MarsdenWeighing](http://www.youtube.com/MarsdenWeighing)

**P** <http://pinterest.com/marsdenweighing>