

MARSDEN

M-900

Marsden Bed Weigher

Easy to use scale for bed-bound patients

Class III Approved and accurate to 200g

Wheels and handles for easy transportation

Includes Tare and BMI functions

www.marsden-weighing.co.uk



Marsden M-900:

An easy to use bed weighing scale for hospitals and care homes

The M-900 is a portable medical scale that makes weighing bed-bound patients easy.

Get a complete picture of a patient's profile

With the M-900, bed-bound patients no longer need to suffer the discomfort of using a floor scale or hoist weigher.

Two easily portable beams can be placed in front of the castors of their bed, and the bed rolls onto the beams. The digital indicator clearly displays the weight, and Tare can be used to deduct the weight of the bed.

For hospitals, care homes and nursing homes moving to telehealth solutions, this bed weigher, like most of our medical scales, has optional Bluetooth connectivity for transferring patient weight data to databases, smartphones, tablets and PCs.

An optional printer is also available.

Suitable for bariatric beds



Optional thermal printer





Easy to use function keys



Carry handles for easy transportation

Easy to use

With the rechargeable battery fitted as standard this scale can be used absolutely anywhere regardless of where the nearest plug socket is. An AC adaptor is also supplied.

The scale has two displays which show weight, height and BMI. By using the full numeric keypad height can be manually entered to calculate the BMI.

The M-900 is available with a 600kg capacity, weighing in 200g increments.

This means the M-900, a medically approved scale, is accurate enough to detect small changes in a patient's weight.

Built to last

Like other scales in the Marsden range the M-900 is built to last. It comes with a two year warranty, and for extra peace of mind this scale can be kept accurate and reliable with a Marsden service contract.

Product Specifications

Capacity: 600kg
Graduations: 200g
Weight of scale: 11.5kg per beam
Power Supply: Mains/Rechargeable Battery
Dimensions: 1240mm x 360mm x 65mm

Optional Extras:

Bluetooth Connectivity
Thermal Printer



Accuracy Assured

Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: sales@marsdengroup.co.uk

www.marsden-weighing.co.uk

 MarsdenWeighing

 @MarsdenWeighing

 www.youtube.com/MarsdenWeighing

 <http://pinterest.com/marsdenweighing>