

MARSDEN

M-200

Marsden High Capacity Chair Scale



Premium, Class III Approved Chair Scale

Designed to accommodate bariatric patients

Robust construction with comfortable one-piece seat

Rechargeable battery: 80 hours' use between charges

www.marsden-weighing.co.uk



Marsden M-200: A high capacity, high quality chair scale

The M-200 is Marsden's flagship chair scale. **Robustly built**, with **premium features** and an **extra-wide seat**, this chair scale will weigh **almost any patient** - providing reliable weighings for many years.

High accuracy, high capacity

The M-200 is available as either a **250kg or 300kg** chair scale - meaning it can be used to weigh almost any patient.

Graduations start at just **50g**, meaning you'll also get a precise reading of your patient's weight. This means this **Class III Approved** scale can be used for medical diagnosis - as well as gaining a weight for critical treatment - for patients unable to stand on a floor scale.



A thought-out design for easy of use

The M-200 boasts a powder-coated frame to make it more **hard-wearing**. The **extra-wide seat** is a one-piece moulding and provides comfort for the patient.

To make the scale easier to get into (and out of), the **armrests and footrests are hinged**. Brakes are provided on two of the wheels.

An **optional handle** is also available for the back of the scale, for improved control when pushing the scale.





Easy to use features

The M-200 features the popular Marsden DP-3810 indicator. This means its features and functionality is identical to that of many of our other medical scales, so it's **quicker to get to grips with**.

A **full numeric keypad** allows you to **calculate BMI** quickly and easily. Tare allows you to remove unwanted weight and Hold stabilises the weight reading on the clear display. Clear **dual displays** mean weight, height and BMI data is easy to read, instantly

Key features

- Class III Approved
- Heavy duty construction
- One-piece extra wide seat
- Hinged armrests/footrests
- Brakes on two wheels
- Large, bright display
- BMI/Hold/Tare

Use it anywhere

The internal rechargeable battery provides 80 hours of use between charges.

For even greater flexibility, optional **Bluetooth** or **wifi** connectivity is available, allowing you to **wirelessly transmit patient data** from the scale to a device or database. An optional **micro USB port** is also available. Please contact us for more information about wireless data transfer.

Specification

- Capacity: 250kg or 300kg
- Graduations: 50g <150kg> 100g
- Power: Mains adaptor or rechargeable battery
- Dimensions: 530mm x 460mm
- Weight of Scale: Approx. 27kg

- Optional Extras
 - Printer (TP-2100)
 - Bluetooth/Wifi/USB
 - Handle
 - Alternative footrests



Marsden M-200: Key features at-a-glance

Easy to use DP-3810 indicator with large displays

Easy access functions: Hold, Tare, Zero, BMI

Highly robust construction

Optional Bluetooth, wifi or USB connectivity

Hinged armrests

Internal rechargeable battery providing 80 hours' use

Brakes on two wheels

One-piece extra wide seat

Hinged footrests



Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: sales@marsdengroup.co.uk

www.marsden-weighing.co.uk

 MarsdenWeighing

 @MarsdenWeighing

 www.youtube.com/MarsdenWeighing

 <http://pinterest.com/marsdenweighing>

MARSDEN