

MARSDEN

H-630 Wall Mounted Height Measure



- Lightweight aluminium construction
- Sliding cursor for accurate measuring
- 1mm increments
- Fits to a wall and slider folds flat when not in use

Marsden H-630: Wall Mounted Height Measure

- Lightweight aluminium construction
- 215cm measuring range
- Sliding cursor for accurate measuring
- 1mm increments
- Fits to a wall and slider folds flat when not in use
- Attaches to a wall with two hinges (supplied)
- 1 year warranty as standard
- Ideal for use in clinics and GP surgeries



Multi-use height measure

The H-630 is a wall mounted height measure which was designed to weigh a broad range of individuals including children, adults and the elderly.

Simply mount the height measure on the wall of your clinic or GP surgery using the two hinges which are supplied as standard, as you are ready to begin getting accurate readings for your patients.

Fold flat cursor

The GP surgery was in the forefront of our mind when we were designing the H-630 height measure and for this reason, the cursor of the height measure folds flat.

This is because we know that space can be tight within surgeries or clinics and to prevent any harm coming to patients or damage to the device, the cursor can be laid flat when not in use.

Easy to use

The M-630 is a simple and easy to use height measure, which takes seconds to use, once it is in position and affixed to the wall where it will be situated.

A patient needs to simply walk up to the height measure, removing shoes if necessary and the healthcare worker will then slowly move the cursor down until it sits on the patients head.

Measuring Range: 75cm to 215cm

Graduations: 1mm

Dimensions: 880mm x 100mm x 330mm

Weight: 0.6kg