


MARSDEN

HM-250P Portable Height Measure




- Lightweight plastic construction
- Includes wall stabilisers
- 1mm increments
- Widely used by health visitors

Marsden HM-250P:

Lightweight Height Measure

- Lightweight design making it ideal for health visitors
- Portable height measure with a carry case as standard
- Includes wall stabilisers
- Measures up to 2m in 1mm increments



Up to 2m

The HM-250P measures up to 2m in height with increments as low as 1mm. It is ideal for use by school nurses, as it can accommodate the height of all children.

The HM-250P is constructed from lightweight plastic, yet it is incredibly durable and will stand the test of time for any health visitor or school nurse which has this height measure in their arsenal.

Extreme accuracy

With 1mm increments, height stopper and red arrows for precision, the Marsden HM-250P is one of the most accurate height measures available on the market.

Whatever the reason the height measurement is needed for, you can guarantee that you will obtain a precise and accurate measurement with the help of the HM-250P.

Lightweight

The HM-250P is an incredibly lightweight construction height measure, which makes it the ideal tool of choice to health visitors and school nurses.

For this reason, we supply the HM-250P as standard with a corrugated plastic case, which enables the user to conveniently transport the carry case thanks to the easy to use handle

Measuring Range: 13.5cm - 210cm

Graduations: 1mm

Dimensions: 2130mm x 350mm x 415mm

Weight: 2kg